

A

www.dk.co.kr

Elpo-32

Portable X-ray System

Elpo-32 is a cost-effective portable X-ray system that guarantees ultimate productivity. It has a designed friendly and compact design for operator.




Features

- High Frequency Portable X-Ray
- User-Friendly Function
- Dual Integrated Laser Pointer
- Light and Compact Design
- Built in Lever for Operator
- Provide high resolution imaging
- Compact size and lightest weight

DK⁺Medical Systems

① Appearance



Description

- ① Handle
- ② Main Body
- ③ Operation Panel
- ④ Collimator
- ⑤ Main Switch
- ⑥ Mono Jack
- ⑦ Power Inlet
- ⑧ Goniometer
- ⑨ Tapeline
- ⑩ Collimator Dial
- ⑪ Hand Switch
- ⑫ Laser Pointer

📄 Specifications

Elpo-32	
Input Voltage	220~240Vac
Max, Output	3.2kW
Kv, mA, mAs range	40~100kV, 25~60mA, 0.32~100mAS
Laser Pointer	Standard

👉 Optional



H.Q / DK Bldg, 18, Baumoe-ro 7-gil, Seocho-gu, Seoul, 06762, Republic of Korea
Tel +82.2.529.6190 Fax +82.2.577.6194

DK⁺Medical Systems Co.,Ltd
www.dk.co.kr / global@dk.co.kr

FACTORY / 52, Chupalsandan 1-gil, Paengseong-eup, Pyeongtaek-si,
Gyeonggi-do, 17998, Republic of Korea
Tel +82.31.658.7231
Fax +82.31.658.7239



B

Elpo-32

Portable X-ray System

Elpo-32 is a cost-effective portable X-ray system that guarantees ultimate productivity. It has a designed friendly and compact design for operator.



“

High Frequency Portable X-Ray

User-Friendly Function

Dual Integrated Laser Pointer

Light and Compact Design

Built in Lever for Operator


Provide high resolution imaging

Compact size and lightest weight

”

DK⁺ Medical Systems

① Appearance



Description

- ① Handle
- ② Main Body
- ③ Operation Panel
- ④ Collimator
- ⑤ Goniometer
- ⑥ Tapeline
- ⑦ Main Switch
- ⑧ Power Inlet
- ⑨ Mono Jack
- ⑩ Collimator Dial
- ⑪ Hand Switch
- ⑫ Laser Pointer

① Specifications

Elpo-32	
Input Voltage	220~240Vac
Max. Output	3.2kW
Kv, mA, mAs range	40~100kV, 25~60mA, 0.32~100mAS
Laser Pointer	Standard

① Optional



DK⁺ Medical Systems
www.dk.co.kr / global@dk.co.kr

H.Q / DK Bldg, 18, Baumoe-ro 7-gil, Seocho-gu, Seoul, 06762, Republic of Korea
Tel +82.2.529.6190 Fax +82.2.577.6194

FACTORY / 52, Chupalsandan 1-gil, Paengseong-eup, Pyeongtaek-si,
Gyeonggi-do, 17998, Republic of Korea
Tel +82.31.658.7231
Fax +82.31.658.7239

